

Digitek – individuelle Abutments, lieferbar für folgende Implantat-Systeme

Marke	System	Plattform	Individuelles Zirkonabutment (Zr/Ti)	Individuelles Titanabutment
Nobe Biocare	Replace	NP	✓	✓
		RP	✓	✓
		WP	✓	✓
		6,0	✓	✓
	Active/Conical	3,0	✓	✓
		NP	✓	✓
		RP	✓	✓
		WP	✗	✓
	Branemark	NP	✓	✓
		RP	✓	✓
		WP	✓	✓
Straumann	Tissue Level	NCC	✓	✓
		RN	✓	✓
		WN	✓	✓
	Bone Level	NC	✓	✓
		RC	✓	✓
Dentsply	AstraTech OsseoSpeed	3,0	✓	✓
		3,5/4,0	✓	✓
		4,5/5,0	✓	✓
	AstraTech EV	3,6	✓	✓
		4,2	✓	✓
		4,8	✓	✓
	Frialit/XIVE	3,4	✓	✓
		3,8	✓	✓
		4,5	✓	✓
		5,5	✓	✓
	Ankylos	C/X	✓	✗
Megagen	Anyridge	–	✓	✓
3i Biomet	Certain	3,4	✓	✓
		4,1	✓	✓
		5,0	✓	✓
	External Hex	3,4	✓	✓
		4,1	✓	✓
		5,0	✓	✓
Zimmer	Tapered Screw-Vent	3,5	✓	✓
		4,5	✓	✓
		5,7	✓	✓
Osstem	ST	Mini	✓	✓
		Standard	✓	✓
Camlog	Camlog	3,8	✓	✓
		4,3	✓	✓
		5,0	✓	✓
AB	AB (AB)	Narrow	✗	✓
		Standard	✗	✓
Medentis	ICX	–	✓	✓
MIS	C1/V3	NP	✓	✓
	C1/V3	SP	✓	✓
	C1	WP	✓	✓
	Seven	NP	✓	✓
		SP	✓	✓
Euroteknika	Nat. (EUN)	Standard	✓	✓
Bego	Semados	3.25/ 3.75	✗	✓
		4,1	✗	✓
		4,5	✗	✓
		5,5	✗	✓